# Medtronic

*Alleviating Pain · Restoring Health · Extending Life*

**Midas Rex® Curved Bur**  
**Fraise courbe Midas Rex®**  
**Midas Rex® Gekrümmter Schleifbohrer**  
**Fresa curvada Midas Rex®**  
**Midas Rex® Curved Bur**  
**Fresa curva Midas Rex®**  
**Kaareva Midas Rex® -poranterä**  
**Midas Rex® svängd borr**  
**Γωνιακή φρέζα Midas Rex®**  
**Midas Rex® kurvebor**  
**Trépano curvo Midas Rex®**  
**Zakrzywione wiertło Midas Rex®**  
**Zakřivená fréza Midas Rex®**  
**Midas Rex® hajlított fúró**  
**Midas Rex® kurvet bor**  
**Midas Rex® Eğri Frez**

Product Information and Instructions · Informations  
et instructions relatives au produit · Informationen und  
Anleitungen zum Produkt · Información del producto  
e instrucciones · Productinformatie en instructies  
· Informazioni sul prodotto e istruzioni · Tuotetiedot  
ja käyttöohjeet · Πληροφορίες προϊόντος και οδηγίες  
· Produktinformation og -vejledning · Informações  
do produto e instruções · Informacje o produkcie  
oraz instrukcje · Údaje o produktu a návod k použití  
· A termék bemutatása és használata · Produktinformasjon  
og instruksjoner · Ürün Bilgisi ve Talimatlar ·

## Device Description

The Midas Rex® Curved Bur is designed to maximize operative site visibility. It is powered by the Legend EHS® and Legend®/MR7 Pneumatic High Speed motors for use with the Legend AT10 (V01) telescoping attachment base.

## Indications

The Midas Rex® Curved Burs are intended for use in surgical procedures for the following medical applications: Neurosurgical; Spine; Ear, Nose, and Throat; Orthopaedic Surgery, and General and Plastic Surgery, including Maxillofacial, Craniofacial, and Sternotomy. The Curved Bur will be used to cut and/or remove bone and biomaterial.

## Contraindications

None

### WARNING

- Always inspect the components for any damage, before and after use. If damage is observed, do not use the part. Damaged parts may deposit debris on the surgical site.
- Do not use the device if the package is opened or damaged.
- Use adequate irrigation and keep the cooling sleeve soaked during dissection. Inadequate irrigation may cause thermal necrosis.
- Do not modify the bur. Bending or prying may break the bur, causing harm to the patient or operating room staff.
- Excessive pressure applied to the bur may cause bur damage. If this occurs, use extreme care to ensure that all fragments of the bur are removed from the patient.
- Disposable devices are for single-use only. Do not attempt to sterilize disposable devices. The disposables are packed sterile and are not intended for repeat use. To prevent contamination, use only once.
- Test for bur wobble (eccentricity) at the desired speed prior to use. Select a new bur or reduce speed if wobble is observed prior to use or during the procedure. Bur wobble may cause patient injury.


### CAUTION

- Curved burs are not designed for variable tool exposure. Do not attempt to adjust tool exposure, as this may damage the device.
- When operating or testing the motor, ensure the bur is properly inserted and locked into AT10 attachment.

## Instructions for Use

### Bur Installation:

1. Turn the tube locking ring and the tool locking ring on the AT10 (V01) Telescoping Attachment into the  position.
2. Slide the curved bur into the AT10 (V01), until the hub is fully seated.
3. Gently press on the head of the dissecting tool to lock it into place. You will hear a click when the curved bur is fully seated.
4. Turn both locking rings to the  position on the AT10 (V01).
5. Gently pull on the curved bur to verify proper installation.
6. Refer to your motor's *Instruction Manual* for information on setting up the AT10 (V01) attachment and motor.
7. Briefly run the motor with the curved bur installed, checking for bur wobble or excessive vibration. If either occurs, perform the corrective actions described in the *Troubleshooting* section of these instructions.

### Cooling the Bur:

1. Prior to initial use, soak the cooling sleeve by dipping it into a cup of saline or DI water, as shown below.

Wetting of the cooling sleeve prior to cutting

2. During use, maintain copious irrigation of the cooling sleeve and bur tip by dribbling saline or DI water along the entire length of the cooling sleeve and bur tip.

### Bur Operation:

1. Activate the motor and gently press the curved bur against the bone to begin dissection.
2. Use a light sweeping motion to continue removing bone.

#### WARNING

Heavy side loads and/or long operating periods may cause the device to overheat. Do not use an overheated device, as it may cause thermal injury to the patient or operator.

### Bur Removal:

1. Turn both AT10 (V01) locking rings to the  position and pull the bur out of the attachment.
2. Discard the curved bur, according to hospital procedures.

## Reuse and Cleaning

Do not reuse.

## Sterility

Each bur is gamma-sterilized and is not intended for repeated use. Do not attempt to re-sterilize.

## Cutting Time

Continuous cutting for extended periods may cause the device to heat to an uncomfortable temperature. To avoid this, limit cutting to the recommended times below:

	Maximum Continuous Cutting Time	Maximum Total Cutting Time
1 & 2 mm burs	2 minutes	12 minutes
3 & 4 mm burs	3 minutes	12 minutes

## Troubleshooting

IF.....Upon inspection, bur appears to be damaged or defective...

Remove and replace defective bur.

IF.....Bur does not rotate or does not cut bone when pressure is applied...

- Make sure the bur is pressed firmly into the drill. You will feel the bur click into place.
- Make sure the tool locking ring is in the  position.
- Make sure the curved bur hub is fully inserted into the AT10 (V01) and the tube locking ring is in the  position.

IF.....Curved Bur has excessive vibration or wobble...

- Inspect the bur head. If it is visibly bent, replace the bur.
- Make sure the tool locking ring is in the  position.

- Make sure the curved bur hub is fully inserted into the AT10(V01) and the tube locking ring is in the  position.

**IF.....Curved bur hub is loose...**

Make sure the hub is fully inserted into the AT10 (V01) and the tube locking ring is in the  position.

**IF.....Bur continues to malfunction...**

- Remove and replace bur.
- Call Customer Service.

## Limited Warranty

- A. This LIMITED WARRANTY provides assurance for the customer who purchases a Midas Rex® Curved Bur (hereinafter the "Product") that should the Product fail to function to Medtronic Powered Surgical Solutions' published specifications prior to the Product expiration date, Medtronic Powered Surgical Solutions will either replace, or issue a credit for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Powered Surgical Solutions or from its affiliate or its authorized distributor or representative.
- B. To qualify for this LIMITED WARRANTY, the following conditions must be met:
  - 1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
  - 2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
  - 3) Medtronic Powered Surgical Solutions must be notified in writing within thirty (30) days following discovery of a defect.
  - 4) The Product must be returned to Medtronic Powered Surgical Solutions within thirty (30) days of Medtronic Powered Surgical Solutions receiving notice as provided for in (3) above.
  - 5) Upon examination of the Product by Medtronic Powered Surgical Solutions, Medtronic Powered Surgical Solutions shall have determined that: (i) the Product was not altered by anyone other than Medtronic Powered Surgical Solutions or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.
- C. This LIMITED WARRANTY is limited to its express terms. **THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** In no event shall Medtronic Powered Surgical Solutions be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected,

and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

\*This Limited Warranty is provided by Medtronic Powered Surgical Solutions, 4620 North Beach Street, Fort Worth, Texas 76137-4116. It applies only in the United States. Areas outside the United States should contact their Medtronic regional distributor for terms of the warranty.

## Returns

Contact Medtronic Powered Surgical Solutions Customer Service at (800-335-9557) to obtain a Return Merchandise Authorization number (RMA#) prior to shipping the Product to Medtronic Powered Surgical Solutions. Please have the original invoice number or purchase order number available to assist in verifying warranty information. The RMA# should be prominently displayed on the box and included on all paperwork enclosed with the return. All Product returned to Medtronic Powered Surgical Solutions should be safely packed in protective wrapping. Customer must supply the Purchase Order number; the correct shipping and billing address; and a statement of the problem or reason for return.

## Customer Service Information

For further information regarding the use of this product or to report any problems, please contact Medtronic Powered Surgical Solutions at (800-335-9557).

## Symbols

Use with

Reference number

Lot number

Attention, See Instructions For Use

Use by date

For single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Sterilized by gamma irradiation.

Rx Only

United States federal law restricts this device to sale by or on the order of a physician.

USA Only

Authorized Representative in the European Community

Manufacturer

Compliant with European Council Directive MDD 93/42/EEC.

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Medtronic B.V.